Gebrauchsanweisung MED DARM REPAIR

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Lagerbedingungen:

Instructions for Use MED DARM REPAIR

Medical device of class IIb
(suitable for long-term use)

Net content:

Powder: 50 g, 100 g, 200 g, 400 g
Capsules: 50 capsules, 100 capsules,
200 capsules

Ingredients: PMA-Zeolite®

Intended Use:

For leaky gut and irritable bowel syndrome.

For the treatment of a limited barrier function of the intestinal wall (leaky gut) and associated inflammatory processes.

MED DARM REPAIR is a natural medical device that has a restoring effect on the gastrointestinal wall and binds harmful substances (lead, cadmium, arsenic, chromium, nickel and ammonium) in the gastrointestinal tract. It is scientifically recognised that different, particularly recurrent, stresses (psychological burden, diseases and drug therapies) can lead to a disturbance of the integrity of the intestinal wall (leaky gut). Possible symptoms associated with leaky gut include gastrointestinal problems, increased liver values, immunodeficiency or fatigue. Please consult your doctor or pharmacist for further information. It has been scientifically proven that PMA-Zeolite® is capable of significantly positively influencing intestinal wall integrity – measurable through a significantly reduced concentration of the 'tight junction modulator' zonulin and the inflammation parameter alpha-1 antitrypsin.

Application:

Powder: Take one teaspoon (1 tsp = 2 g) 1 – 3 times daily mixed in a glass of water (at least 100 ml) with meals. For acute diarrhea, take 3 teaspoons (= 6 g). The maximum daily intake should not exceed 9 g.

Capsules: Take 3 capsules 2 – 3 times a day with a glass of water (at least 100 ml) with meals. For acute diarrhoea, take up to 15 capsules.

The duration of intake should be at least 12 weeks or as long as the recurrent stresses exist. The intake is recommended for adults from the age of 18 up.

Warnings and safety instructions:

There are no undesirable side effects apart of occasional constipation. This is particularly important in case of chronic constipation. Discontinue use if you notice an unexpected effect and consult your physician. As there are not yet sufficient data available, its use cannot be recommended for children and adolescents (under 18 years of age), pregnant and lactating women. To avoid adverse drug interactions, we recommend taking MED DARM REPAIR and medication with a 2-hour time lag, whereby medication should be taken before the PMA-Zeolite®, unless otherwise prescribed. Do not inhale the powder. In the event of serious incidents, please contact the manufacturer and the competent authority. Keep out of reach of children under 6 years of age. Firmly close can after use. Dispose of the device with household waste in accordance with local regulations. Constant quality control ensures high quality and effectiveness as an approved medical device in acc. with European Regulation 2017/745

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PANACEO International GmbH